

ApoC-III Polyclonal Antibody

Catalog No: #41916

Package Size: #41916-1 50ul #41916-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ApoC-III Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ApoC-III Polyclonal Antibody detects endogenous levels of ApoC-III protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human ApoC-III.
Target Name	ApoC-III
Other Names	APOC3; Apolipoprotein C-III; Apo-CIII; ApoC-III; Apolipoprotein C3
Accession No.	Swiss-Prot: P02656NCBI Gene ID: 345
Uniprot	P02656
GeneID	345;
SDS-PAGE MW	11kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T-UV cells using ApoC-III Polyclonal Antibody

Note: This product is for in vitro research use only