

FGF-4 Polyclonal Antibody

Catalog No: #41918

Package Size: #41918-1 50ul #41918-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGF-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	FGF-4 Polyclonal Antibody detects endogenous levels of FGF-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-4.
Target Name	FGF-4
Other Names	FGF4; HST; HSTF1; KS3; Fibroblast growth factor 4; FGF-4; Heparin secretory-transforming protein 1; HST; HST-1; HSTF-1; Heparin-binding growth factor 4; HBGF-4; Transforming protein KS3
Accession No.	Swiss-Prot: P08620NCBI Gene ID: 2249
Uniprot	P08620
GeneID	2249;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 HeLa cells using FGF-4 Polyclonal Antibody

Note: This product is for in vitro research use only