

Stra8 Polyclonal Antibody

Catalog No: #41923

Package Size: #41923-1 50ul #41923-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stra8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Stra8 Polyclonal Antibody detects endogenous levels of Stra8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Stra8.
Target Name	Stra8
Other Names	STRA8; Stimulated by retinoic acid gene 8 protein homolog
Accession No.	Swiss-Prot: Q7Z7C7NCBI Gene ID: 346673
Uniprot	Q7Z7C7
GeneID	346673;
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3, 293 cells using Stra8 Polyclonal Antibody.

Note: This product is for in vitro research use only