

RGS1 Polyclonal Antibody

Catalog No: #41925

Package Size: #41925-1 50ul #41925-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RGS1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	RGS1 Polyclonal Antibody detects endogenous levels of RGS1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human RGS1.
Target Name	RGS1
Other Names	RGS1; 1R20; BL34; IER1; Regulator of G-protein signaling 1; RGS1; B-cell activation protein BL34; Early response protein 1R20
Accession No.	Swiss-Prot: Q08116NCBI Gene ID: 5996
Uniprot	Q08116
GeneID	5996;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using RGS1 Polyclonal Antibody.

Note: This product is for in vitro research use only