

GHR Polyclonal Antibody

Catalog No: #41939

Package Size: #41939-1 50ul #41939-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GHR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	GHR Polyclonal Antibody detects endogenous levels of GHR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GHR.
Target Name	GHR
Other Names	GHR; Growth hormone receptor; GH receptor; Somatotropin receptor
Accession No.	Swiss-Prot: P10912NCBI Gene ID: 2690
Uniprot	P10912
GeneID	2690;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of SKOV3 cells using GHR Polyclonal Antibody.

Note: This product is for in vitro research use only