

CD96 Polyclonal Antibody

Catalog No: #41963

Package Size: #41963-1 50ul #41963-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD96 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD96 Polyclonal Antibody detects endogenous levels of CD96 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD96.
Target Name	CD96
Other Names	CD96; T-cell surface protein tactile; Cell surface antigen CD96; T cell-activated increased late expression protein; CD96
Accession No.	Swiss-Prot: P40200NCBI Gene ID: 10225
Uniprot	P40200
GeneID	10225;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of 293-UV cells using CD96 Polyclonal Antibody

Western Blot analysis of 293-UV cells using CD96 Polyclonal Antibody



Note: This product is for in vitro research use only