

## Nectin 1 Polyclonal Antibody

Catalog No: #41967

Package Size: #41967-1 50ul #41967-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nectin 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Nectin 1 Polyclonal Antibody detects endogenous levels of Nectin 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Nectin 1.
Target Name	Nectin 1
Other Names	PVRL1; HVEC; PRR1; Poliovirus receptor-related protein 1; Herpes virus entry mediator C; Herpesvirus entry mediator C; HveC; Herpesvirus Ig-like receptor; HlgR; Nectin-1; CD111
Accession No.	Swiss-Prot: Q15223NCBI Gene ID: 5818
Uniprot	Q15223
GeneID	5818;
SDS-PAGE MW	57kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of KB cells using Nectin 1 Polyclonal Antibody



Western Blot analysis of KB cells using Nectin 1 Polyclonal Antibody

Note: This product is for in vitro research use only