

IL-9R Polyclonal Antibody

Catalog No: #41971

Package Size: #41971-1 50ul #41971-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-9R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-9R Polyclonal Antibody detects endogenous levels of IL-9R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-9R.
Target Name	IL-9R
Other Names	IL9R; Interleukin-9 receptor; IL-9 receptor; IL-9R; CD129
Accession No.	Swiss-Prot: Q01113NCBI Gene ID: 3581
Uniprot	Q01113
GeneID	3581;
SDS-PAGE MW	58kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of KB cells using IL-9R Polyclonal Antibody



Western Blot analysis of KB NIH-3T3 cells using IL-9R Polyclonal Antibody

Note: This product is for in vitro research use only