

Ox40 Polyclonal Antibody

Catalog No: #41972

Package Size: #41972-1 50ul #41972-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ox40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Ox40 Polyclonal Antibody detects endogenous levels of Ox40 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ox40.
Target Name	Ox40
Other Names	TNFRSF4; TXGP1L; Tumor necrosis factor receptor superfamily member 4; ACT35 antigen; OX40L receptor; TAX transcriptionally-activated glycoprotein 1 receptor; CD134
Accession No.	Swiss-Prot: P43489NCBI Gene ID: 7293
Uniprot	P43489
GeneID	7293;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of PC-12 cells using Ox40 Polyclonal Antibody

Note: This product is for in vitro research use only