

Blr1 Polyclonal Antibody

Catalog No: #41978

Package Size: #41978-1 50ul #41978-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Blr1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Blr1 Polyclonal Antibody detects endogenous levels of Blr1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Blr1.
Target Name	Blr1
Other Names	CXCR5; BLR1; MDR15; C-X-C chemokine receptor type 5; CXC-R5; CXCR-5; Burkitt lymphoma receptor 1; Monocyte-derived receptor 15; MDR-15; CD185
Accession No.	Swiss-Prot: P32302NCBI Gene ID: 643
Uniprot	P32302
GeneID	643;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using Blr1 Polyclonal Antibody

Western Blot analysis of NIH-3T3 cells using Blr1 Polyclonal Antibody



Note: This product is for in vitro research use only