

IL-12R β 1 Polyclonal Antibody

Catalog No: #41979

Package Size: #41979-1 50ul #41979-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IL-12R β 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-12R β 1 Polyclonal Antibody detects endogenous levels of IL-12R β 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-12R β 1.
Target Name	IL-12R β 1
Other Names	IL12RB1; IL12R; IL12RB; Interleukin-12 receptor subunit beta-1; IL-12 receptor subunit beta-1; IL-12R subunit beta-1; IL-12R-beta-1; IL-12RB1; IL-12 receptor beta component; CD212
Accession No.	Swiss-Prot: P42701NCBI Gene ID: 3594
Uniprot	P42701
GeneID	3594;
SDS-PAGE MW	73kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of HeLa cells using IL-12R β 1 Polyclonal Antibody

Note: This product is for in vitro research use only