

EMR2 Polyclonal Antibody

Catalog No: #41989

Package Size: #41989-1 50ul #41989-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EMR2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	EMR2 Polyclonal Antibody detects endogenous levels of EMR2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human EMR2.
Target Name	EMR2
Other Names	EMR2; EGF-like module-containing mucin-like hormone receptor-like 2; EGF-like module receptor 2; CD312
Accession No.	Swiss-Prot: Q9UHX3NCBI Gene ID: 30817
Uniprot	Q9UHX3
GeneID	30817;
SDS-PAGE MW	72kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF-7 cells using EMR2 Polyclonal Antibody

Note: This product is for in vitro research use only