

IL-10R α Polyclonal Antibody

Catalog No: #41994

Package Size: #41994-1 50ul #41994-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IL-10R α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-10R α Polyclonal Antibody detects endogenous levels of IL-10R α protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-10R α .
Target Name	IL-10R α
Other Names	IL10RA; IL10R; Interleukin-10 receptor subunit alpha; IL-10 receptor subunit alpha; IL-10R subunit alpha; IL-10RA; CDw210a; Interleukin-10 receptor subunit 1; IL-10R subunit 1; IL-10R1; CD210
Accession No.	Swiss-Prot: Q13651NCBI Gene ID: 3587
Uniprot	Q13651
GeneID	3587;
SDS-PAGE MW	63kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of PC12 cells using IL-10R α Polyclonal Antibody

Note: This product is for in vitro research use only