

Amphiregulin Polyclonal Antibody

Catalog No: #42008

Package Size: #42008-1 50ul #42008-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amphiregulin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Amphiregulin Polyclonal Antibody detects endogenous levels of Amphiregulin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Amphiregulin.
Target Name	Amphiregulin
Other Names	AREG; SDGF; AREGB; Amphiregulin; AR; Colorectum cell-derived growth factor; CRDGF
Accession No.	Swiss-Prot: P15514NCBI Gene ID: 374/727738
Uniprot	P15514
GeneID	374;
SDS-PAGE MW	29kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HBE cells using Amphiregulin Polyclonal Antibody

Note: This product is for in vitro research use only