

Protein AMBP Monoclonal Antibody

Catalog No: #42043

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Protein AMBP Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	protein G purified
Applications	IHC
Species Reactivity	Hu
Specificity	specific for Human protein AMBP denatured and native forms
Immunogen Type	protein
Immunogen Description	Recombinant Human protein AMBP
Target Name	Protein AMBP
Other Names	Protein HC, Alpha-1 microglycoprotein, ITI-LC, Bikunin, EDC1, HI-30
Accession No.	Swiss-Prot#: P02760
Uniprot	P02760
GeneID	259;
Concentration	1.0mg/mL
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human liver tissue using #42043 at dilution of 1:200.

Background

Inter-alpha-trypsin inhibitor inhibits trypsin, plasmin, and lysosomal granulocytic elastase. Inhibits calcium oxalate crystallization. Trypstatin is a trypsin inhibitor

References

- [1] "Sequence of a full length cDNA coding for human protein HC (alpha 1 microglobulin)." Traboni C., Cortese R. Nucleic Acids Res. 14:6340-6340(1986) [2] "The mRNA for a proteinase inhibitor related to the HI-30 domain of inter-alpha-tr

Note: This product is for in vitro research use only