

# Human Tumor necrosis factor receptor superfamily member 3 (LTBR) ELISA Kit



Catalog No: #EK10014

Orders: order@signalwayantibody.com

Package Size: #EK10014-1 48T #EK10014-2 96T

Support: tech@signalwayantibody.com

## Description

Product Name	Human Tumor necrosis factor receptor superfamily member 3 (LTBR) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	D12S370; LT-BETA-R; TNF-R-III; TNFCR; TNFR-RP; TNFR2-RP; TNFRSF3; TNFR superfamily; member 3 lymphotoxin B receptor lymphotoxin beta receptor tumor necrosis factor C receptor tumor necrosis fa
Accession No.	P36941
Uniprot	P36941
GeneID	4055;
Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).

## Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.078 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 &amp;mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:Sandwich

**Test principle:**This assay employs a two-site sandwich ELISA to quantitate LTBR in samples. An antibody specific for LTBR has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyLTBR present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for LTBR is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of LTBR bound in the initial step. The color development is stopped and the intensity of the color is measured.

**Product Overview:**This protein is a member of the tumor necrosis factor (TNF) family of receptors. It is expressed on the surface of most cell types, including cells of epithelial and myeloid lineages, but not on T and B lymphocytes. The protein specifically binds the lymphotoxin membrane form (a complex of lymphotoxin-alpha and lymphotoxin-beta). The encoded protein and its ligand play a role in the development and organization of lymphoid tissue and transformed cells. Activation of the encoded protein can trigger apoptosis.

---

Note: This product is for in vitro research use only