

## human Protein AMBP Polyclonal Antibody

Catalog No: #42062

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	human Protein AMBP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total human Protein AMBP polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Protein AMBP protein
Target Name	human Protein AMBP
Other Names	Alternative name(s): Alpha-1 microglycoprotein Complex-forming glycoprotein heterogeneous in charge Bikunin EDC1 HI-30 Uronic-acid-rich protein AMBP HCP, ITIL
Accession No.	Swiss-Prot#: P02760
Uniprot	P02760
GeneID	259;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human placenta using #42062 at dilution of 1:20.

## Background

Inter-alpha-trypsin inhibitor inhibits trypsin, plasmin, and lysosomal granulocytic elastase. Inhibits calcium oxalate crystallization.

## References

---

[1]Sequence of a full length cDNA coding for human protein HC (alpha 1 microglobulin).Traboni C., Cortese R.Nucleic Acids Res. 14:6340-6340(1986)

[2]The mRNA for a proteinase inhibitor related to the HI-30 domain of inter-alpha-trypsin inhibitor also en

---

Note: This product is for in vitro research use only