

α 1-acidglycoprotein Polyclonal Antibody

Catalog No: #42087

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	α 1-acidglycoprotein Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total α 1-acidglycoprotein polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Human α 1-acidglycoprotein protein
Target Name	α 1-acidglycoprotein
Other Names	α 1-AGP
Calculated MW	45kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: 1:500 - 1:1000

Images



All lanes: α 1-acidglycoprotein antibody at 2ug/ml+Human serum at 1 : 50 secondary
Goat polyclonal to rabbit at 1/10000 dilution
predicted band size :45kDa
observed band size :30,50kDa

Note: This product is for in vitro research use only