

## CDKN2AIP N-terminal-like protein Polyclonal Antibody

Catalog No: #42115

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CDKN2AIP N-terminal-like protein Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total CDKN2AIP N-terminal-like protein polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Human CDKN2AIP N-terminal-like protein protein
Target Name	CDKN2AIP N-terminal-like protein
Other Names	CDKN2A-interacting protein N-terminal-like protein, CDKN2AIPNL
Accession No.	Swiss-Prot#: Q96HQ2
Uniprot	Q96HQ2
GeneID	91368;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded mouse liver using #42115 at dilution of 1:100.



Immunohistochemical analysis of paraffin-embedded colon cancer using #42115 at dilution of 1:100.

## References

[1] Mural R.J., Istrail S., Sutton G.G., Florea L., Halpern A.L., Mobarry C.M., Lippert R., Walenz B., Shatkay H., Dew I., Miller J.R., Flanigan M.J., Edwards N.J., Bolanos R., Fasulo D., Halldorsson B.V., Hannenhalli S., Turner R. , Yooseph S., Lu F., N

Note: This product is for in vitro research use only