

D-amino acid oxidase activator Polyclonal Antibody

Catalog No: #42143

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	D-amino acid oxidase activator Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total D-amino acid oxidase activator polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human D-amino acid oxidase activator protein
Target Name	D-amino acid oxidase activator
Other Names	Protein G72 DAOA G72
Accession No.	Swiss-Prot#: P59103
Uniprot	P59103
GeneID	267012;
Calculated MW	18kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes: D-amino acid oxidase activator antibody at 2ug/ml

Lane 1: PC-3 whole cell lysate

Lane 2: mouse gonad tissue

Secondary

Goat polyclonal to Rabbit IgG at 1/10000 dilution

Predicted band size: 18kDa

Observed band size: 18kDa

Background

Seems to activate D-amino acid oxidase.

References

[1]Genetic and physiological data implicating the new human gene G72 and the gene for D-amino acid oxidase in schizophrenia.Chumakov I., Blumenfeld M., Guerassimenko O., Cavarec L., Palicio M., Abderrahim H., Bougueleret L., Barry C., Tanaka H., La Rosa P

Note: This product is for in vitro research use only