

Protein delta homolog 1 Polyclonal Antibody

Catalog No: #42148

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Protein delta homolog 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Protein delta homolog 1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Protein delta homolog 1 protein
Target Name	Protein delta homolog 1
Other Names	pG2, Fetal antigen 1, DLK1, DLK
Accession No.	Swiss-Prot#: P80370
Uniprot	P80370
GeneID	8788;
Calculated MW	41kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes: protein delta homolog 1 antibody at 2ug/ml+U251
whole cell lysate
Secondary
Goat polyclonal to Rabbit IgG at 1/10000 dilution
Predicted band size:41kDa
Observed band size:41kDa

Background

May have a role in neuroendocrine differentiation.

References

[1]dlk, a putative mammalian homeotic gene differentially expressed in small cell lung carcinoma and neuroendocrine tumor cell line.Laborda J., Sausville E.A., Hoffman T., Notario V.J. Biol. Chem. 268:3817-3820(1993)[2]dlk, pG2 and Pref-1 mRNAs encode s

Note: This product is for in vitro research use only