

## Annexin A6 Polyclonal Antibody

Catalog No: #42163

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Annexin A6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Annexin A6 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Human Annexin A6 protein
Target Name	Annexin A6
Other Names	67 kDa calelectrin, Annexin VI, Annexin-6, Calphobindin-II, Chromobindin-20, Lipocortin VI, Protein III, p68, p70, ANXA6 , ANX6
Accession No.	Swiss-Prot#: P08133
Uniprot	P08133
GeneID	309;
Calculated MW	76kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500 - 1:1000

## Images



All lanes: Annexin A6 antibody at 2ug/ml+Rat liver tissue  
 Secondary  
 Goat polyclonal to Rabbit IgG at 1/10000 dilution  
 Predicted band size:76kDa  
 Observed band size:76kDa

## Background

May associate with CD21. May regulate the release of Ca<sup>2+</sup> from intracellular stores.

## References

---

[1] Primary structure of the human, membrane-associated Ca<sup>2+</sup>-binding protein p68 a novel member of a protein family. Crompton M.R., Owens R.J., Totty N.F., Moss S.E., Waterfield M.D., Crompton M.J. EMBO J. 7:21-27(1988) [2] Human 67-kDa callectrin contains

---

Note: This product is for in vitro research use only