

# Glycosyltransferase 8 domain-containing protein 2 Polyclonal Antibody

Catalog No: #42185

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Glycosyltransferase 8 domain-containing protein 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total Glycosyltransferase 8 domain-containing protein 2 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Glycosyltransferase 8 domain-containing protein 2 protein
Target Name	Glycosyltransferase 8 domain-containing protein 2
Other Names	GLT8D2 GALA4A UNQ1901/PRO4347
Accession No.	Swiss-Prot#: Q9H1C3
Uniprot	Q9H1C3
GeneID	83468;
Calculated MW	40kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting: □1:500 - 1:1000

## Images



All lanes: Glycosyltransferase 8 domain-containing protein 2 antibody at 2ug/ml  
 Lane 1: mouse heart tissue  
 Lane 2: mouse kidney tissue  
 secondary  
 Goat polyclonal to rabbit at 1/10000 dilution  
 predicted band size :40kDa  
 observed band size :40kDa

## References

[1]Gala4a, glycosyltransferase. Steffensen R., Bennett E.P. Submitted (JAN-2000) [2]The secreted protein discovery initiative (SPDI), a large-scale

---

Note: This product is for in vitro research use only