

Mouse Loricrin (LOR) ELISA Kit

Catalog No: #EK11524



Package Size: #EK11524-1 48T #EK11524-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Loricrin (LOR) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse (<i>Mus musculus</i>)
Other Names	MGC111513;
Accession No.	P18165
Uniprot	P18165
GeneID	16939;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:7.8-500 ng/mL

Sensitivity:3.2 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate LOR in samples. An antibody specific for LOR has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyLOR present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for LOR is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of LOR bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**Loricrin is a major protein component of the cornified cell envelope found in terminally differentiated epidermal cells. Mutations in this gene are associated with Vohwinkel's syndrome and progressive symmetric erythrokeratoderma, both inherited skin diseases. There were multiple sequence variants within these 2 size class alleles due to various deletions of 12 bp (4 amino acids) in the major loop of this glycine loop domain. By use of a specific loricrin antibody in immunogold electron microscopy, they showed that loricrin appears initially in the granular layer of human epidermis and forms composite keratohyalin granules with profilaggrin, but localizes to the cell periphery (cell envelope) of fully differentiated stratum corneum cells.

Note: This product is for in vitro research use only