

Heterogeneous nuclear ribonucleoprotein A1 Polyclonal Antibody

Catalog No: #42208

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Heterogeneous nuclear ribonucleoprotein A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Heterogeneous nuclear ribonucleoprotein A1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Heterogeneous nuclear ribonucleoprotein A1 protein
Target Name	Heterogeneous nuclear ribonucleoprotein A1
Other Names	Helix-destabilizing protein, Single-strand RNA-binding protein, hnRNP core protein A1, HNRNPA1, HNRPA1
Accession No.	Swiss-Prot#: P09651
Uniprot	P09651
GeneID	3178;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded mouse heart using #42208 at dilution of 1:100.

Background

Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection. May play a role in HCV RNA replication.

References

[1] Isolation of an active gene encoding human hnRNP protein A1. Evidence for alternative splicing. Biamonti G., Buvoli M., Bassi M.T., Morandi C., Cobianchi F., Riva S.J. Mol. Biol. 207:491-503(1989) [2] cDNA cloning of human hnRNP protein A1 reveals t

Note: This product is for in vitro research use only