

Human Muscarinic acetylcholine receptor M3 (CHRM3) ELISA Kit

Catalog No: #EK11641

Orders: order@signalwayantibody.com

Package Size: #EK11641-1 48T #EK11641-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Muscarinic acetylcholine receptor M3 (CHRM3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	HM3; m3 muscarinic receptor muscarinic acetylcholine receptor M3
Accession No.	P20309
Uniprot	P20309
GeneID	1131;
Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).

Application Details

Detect Range:0.31 ng/mL - 20 ng/mL.

Sensitivity:0.156 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:Sandwich

Test principle:This assay employs a two-site sandwich ELISA to quantitate CHRM3 in samples. An antibody specific for CHRM3 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyCHRM3 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for CHRM3 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of CHRM3 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Product Overview:The , also known as the cholinergic receptor, muscarinic 3, is a muscarinic acetylcholine receptor.The M4 muscarinic receptors are located at many places in the body, e.g. smooth muscles, endocrine, in exocrine glands, as well as in the lungs. They are al????

Note: This product is for in vitro research use only