

39S ribosomal protein L19, mitochondrial Polyclonal Antibody

Catalog No: #42259

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	39S ribosomal protein L19, mitochondrial Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 39S ribosomal protein L19, mitochondrial polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 39S ribosomal protein L19, mitochondrial protein
Target Name	39S ribosomal protein L19, mitochondrial
Other Names	39S ribosomal protein L15, mitochondrial MRPL19 KIAA0104, MRPL15
Accession No.	Swiss-Prot#: P49406
Uniprot	P49406
GeneID	9801;
Calculated MW	34kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details


Western blotting: □ 1:500 - 1:1000

Immunohistochemistry: 1:20 - 1:200

Images



All lanes: 39S ribosomal protein L19, mitochondrial antibody at 2ug/ml
Lane 1: MCF-7 whole cell lysate
Lane 2: mouse liver tissue
Secondary
Goat polyclonal to Rabbit IgG at 1/10000 dilution
Predicted band size: 34kDa
Observed band size: 34kDa



Immunohistochemical analysis of paraffin-embedded human placenta using #42259 at dilution of 1:50.

References

[1]Generation and annotation of the DNA sequences of human chromosomes 2 and 4.Hillier L.W., Graves T.A., Fulton R.S., Fulton L.A., Pepin K.H., Minx P., Wagner-McPherson C., Layman D., Wylie K., Sekhon M., Becker M.C., Fewell G.A., Delehaunty K.D., Miner

Note: This product is for in vitro research use only