

General Abscisic Acid (ABA) ELISA Kit

Catalog No: #EK12370



Package Size: #EK12370-1 48T #EK12370-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	General Abscisic Acid (ABA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	General
Other Names	AB-A; abscisin II; Dormin
Storage	12 months

Application Details

Detect Range:1.56-100ng/mL

Sensitivity:0.77ng/mL

Sample Type:serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only