

14 kDa phosphohistidine phosphatase Polyclonal Antibody

Catalog No: #42291

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	14 kDa phosphohistidine phosphatase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 14 kDa phosphohistidine phosphatase polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 14 kDa phosphohistidine phosphatase protein
Target Name	14 kDa phosphohistidine phosphatase
Other Names	Phosphohistidine phosphatase 1 Protein janus-A homolog PHPT1 PHP14 CGI-202, HSPC141
Accession No.	Swiss-Prot#: Q9NRX4
Uniprot	Q9NRX4
GeneID	29085;
Calculated MW	13kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes: 14 kDa phosphohistidine phosphatase antibody at
 2ug/ml+MCF-7 whole cell lysate
 secondary
 Goat polyclonal to rabbit at 1/10000 dilution
 predicted band size :13kDa
 observed band size :13kDa

Background

Exhibits phosphohistidine phosphatase activity.

References

[1]Identification and characterization of a mammalian 14-kDa phosphohistidine phosphatase.Ek P., Pettersson G., Ek B., Gong F., Li J.-P., Zetterqvist O.Eur. J. Biochem. 269:5016-5023(2002) [2]Gene expression profiling in the human hypothalamus-pituitar

Note: This product is for in vitro research use only