

Putative phospholipase B-like 2 Polyclonal Antibody

Catalog No: #42295

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Putative phospholipase B-like 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Ms
Specificity	The antibody detects endogenous level of total Putative phospholipase B-like 2 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant mouse Putative phospholipase B-like 2 protein
Target Name	Putative phospholipase B-like 2
Other Names	66.3 kDa protein 76 kDa protein LAMA-like protein 2 Lamina ancestor homolog 2 Phospholipase B domain-containing protein 2 Plbd2
Accession No.	Swiss-Prot#: Q3TCN2
Uniprot	Q3TCN2
GeneID	71772;
Calculated MW	66kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes: Putative phospholipase B-like 2 antibody at 2ug/ml+mouse brain tissue
Secondary
Goat polyclonal to Rabbit IgG at 1/15000 dilution
Predicted band size: 66KDa
Observed band size: 66KDa

Background

Putative phospholipase

References

[1]"Positional cloning of the murine flavivirus resistance gene."Perelygin A.A., Scherbik S.V., Zhulin I.B., Stockman B.M., Li Y., Brinton M.A. Proc. Natl. Acad. Sci. U.S.A. 99:9322-9327(2002) [2]"The transcriptional landscape of the

Note: This product is for in vitro research use only