

ES1 protein homolog, mitochondrial Polyclonal Antibody

Catalog No: #42296

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ES1 protein homolog, mitochondrial Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total ES1 protein homolog, mitochondrial polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human ES1 protein homolog, mitochondrial protein
Target Name	ES1 protein homolog, mitochondrial
Other Names	Protein GT335 Protein KNP-I
Accession No.	Swiss-Prot:P30042Gene ID:8209
Uniprot	A0A0B4J2D5
GeneID	8209;
Calculated MW	29kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes : ES1 protein homolog, mitochondrial antibody at at 2 ug/ml
 Lane 1 : EC109 whole cell lysate
 Lane 2 : 293T whole cell lysate
 Secondary
 Goat polyclonal to Rabbit IgG at 1/15000 dilution
 Predicted band size :29kDa
 Observed band size:29kDa

References

[1] "Isolation of cDNA for a novel human protein KNP-I that is homologous to the E. coli SCRP-27A protein from the autoimmune polyglandular disease type I (APECED) region of chromosome 21q22.3." Nagamine K., Kudoh J., Minoshima S., Kawasaki K., Asakawa S

Note: This product is for in vitro research use only