

Monkey Alpha-Fetoprotein (aFP) ELISA Kit

Catalog No: #EK12500



Package Size: #EK12500-1 48T #EK12500-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Monkey Alpha-Fetoprotein (aFP) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rhesus monkey Simian
Other Names	A-FP; FETA; HPAFP; Alpha-Fetoglobulin; Alpha-1-fetoprotein
Storage	12 months

Application Details

Detect Range:1.56-100ng/mL

Sensitivity:1.2ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only