

Chicken Angiotensin II (AngII) ELISA Kit

Catalog No: #EK12591



Package Size: #EK12591-1 48T #EK12591-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken Angiotensin II (AngII) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken Gallus
Other Names	Ang-II; Angiotensin-2
Storage	12 months

Application Details

Detect Range: 31.2-2000pg/mL

Sensitivity: 12.1pg/mL

Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only