

Equine Angiotensin II (AngII) ELISA Kit

Catalog No: #EK12592



Package Size: #EK12592-1 48T #EK12592-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Equine Angiotensin II (AngII) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Equus caballus; Equine Horse
Other Names	Ang-II; Angiotensin-2
Storage	12 months

Application Details

Detect Range: 31.2-2000pg/mL

Sensitivity: 13.3pg/mL

Sample Type: serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 100 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only