

60S ribosomal protein L30 Polyclonal Antibody

Catalog No: #42313

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	60S ribosomal protein L30 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 60S ribosomal protein L30 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 60S ribosomal protein L30 protein
Target Name	60S ribosomal protein L30
Other Names	RPL30
Accession No.	Swiss-Prot#: P62888
Uniprot	P62888
GeneID	6156;
Calculated MW	13kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: 1:500 - 1:1000

Images



All lanes: 60S ribosomal protein L30 antibody at 2ug/ml+A549 whole cell lysate
 Secondary
 Goat polyclonal to Rabbit IgG at 1/15000 dilution
 Predicted band size:13kda
 Observed band size:13kda

References

[1]"Crystallization and preliminary X-ray structure analysis of human ribosomal protein L30e."Kawaguchi A., Ose T., Yao M., Tanaka I.Acta Crystallogr. F 67:1516-1518(2011).[2]"N-terminal acetylome analyses and functional insights of th

Note: This product is for in vitro research use only