

Non-specific lipid-transfer protein Polyclonal Antibody

Catalog No: #42318

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Non-specific lipid-transfer protein Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Non-specific lipid-transfer protein polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Non-specific lipid-transfer protein
Target Name	Non-specific lipid-transfer protein
Other Names	Propanoyl-CoA C-acyltransferase SCP-chi SCPX Sterol carrier protein 2 Sterol carrier protein X SCP2
Accession No.	Swiss-Prot#: P22307
Uniprot	P22307
GeneID	6342;
Calculated MW	59kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes: Non-specific lipid-transfer protein antibody at 2ug/ml
 Lane 1: A549 whole cell lysate
 Lane 2: MCF-7 whole cell lysate
 secondary
 Goat polyclonal to rabbit at 1/10000 dilution
 predicted band size :59kDa
 observed band size :59kDa

Background

Mediates in vitro the transfer of all common phospholipids, cholesterol and gangliosides between membranes. May play a role in regulating steroidogenesis.

References

[1]The structure of the human sterol carrier protein X/sterol carrier protein 2 gene (SCP2).Ohba T., Rennert H., Pfeifer S.M., He Z., Yamamoto R., Holt J.A., Billheimer J.T., Strauss J.F. III Genomics 24:370-374(1994) [2]cDNAs encoding members of a fam

Note: This product is for in vitro research use only