

Transcription initiation protein SPT3 homolog Polyclonal Antibody

Catalog No: #42337

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transcription initiation protein SPT3 homolog Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Transcription initiation protein SPT3 homolog polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Transcription initiation protein SPT3 homolog protein
Target Name	Transcription initiation protein SPT3 homolog
Other Names	SPT3-like protein SUPT3H SPT3
Accession No.	Swiss-Prot#: O75486
Uniprot	O75486
GeneID	8464;
Calculated MW	40kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes : Transcription initiation protein SPT3 homolog at 2ug/ml

Lane 1:mouse brain tissue

Lane 2:colo320 whole cell lysate

Secondary

Goat polyclonal to Rabbit IgG at 1/10000 dilution

Predicted band size:40kDa

Observed band size:40kDa

Background

Probable transcriptional activator.

References

[1]Characterization of a human homologue of the *Saccharomyces cerevisiae* transcription factor Spt3 (SUPT3H). Yu J., Madison J.M., Mundlos S., Winston F., Olsen B.R. *Genomics* 53:90-96(1998)[2]Histone-like TAFs within the PCAF histone acetylase complex. Ogrý

Note: This product is for in vitro research use only