

Human Complement C3 Convertase (C3c) ELISA Kit

Catalog No: #EK12941



Package Size: #EK12941-1 48T #EK12941-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Complement C3 Convertase (C3c) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Homo sapiens human
Other Names	C4b2a
Storage	12 months

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.37ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only