

Human Complement Component 4d (C4d) ELISA Kit

Catalog No: #EK12954



Package Size: #EK12954-1 48T #EK12954-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Complement Component 4d (C4d) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Homo sapiens human
Storage	12 months

Application Details

Detect Range:9.37-600pg/mL

Sensitivity:4.2pg/mL

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only