

Human Complement Component 4 (C4) ELISA Kit

Catalog No: #EK12955



Package Size: #EK12955-1 48T #EK12955-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Complement Component 4 (C4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Homo sapiens human
Other Names	C4BP-a; Complement Component 4 Binding Protein,Alpha
Storage	12 months

Application Details

Detect Range:4.68-300ng/mL

Sensitivity:1.9ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only