

Rat Complement Component 4 (C4) ELISA Kit

Catalog No: #EK12957



Package Size: #EK12957-1 48T #EK12957-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Rat Complement Component 4 (C4) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rattus norvegicus Rat |
| Other Names | C4BP-a; Complement Component 4 Binding Protein,Alpha |
| Accession No. | UniProtKB:Q63514Gene ID:24235 |
| Uniprot | Q63514 |
| GeneID | 24235; |
| Storage | 12 months |

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:25.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only