

## 60S ribosomal protein L15 Polyclonal Antibody

Catalog No: #42354

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	60S ribosomal protein L15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 60S ribosomal protein L15 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 60S ribosomal protein L15 protein
Target Name	60S ribosomal protein L15
Other Names	RPL15
Accession No.	Swiss-Prot#: P61313
Uniprot	P61313
GeneID	6138;
Calculated MW	22.4kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting: □1:500 - 1:1000

## Images



All lanes : 60S ribosomal protein L15 antibody at at 2ug/ml  
 Lane 1 : EC109 whole cell lysate  
 Lane 2 : 293T whole cell lysate  
 Secondary Goat polyclonal to Rabbit IgG at 1/15000 dilution  
 Predicted band size : 22.4kDa  
 Observed band size: 70 kDa

## Background

RPL15 belongs to the ribosomal protein L15e family.

## References

---

[1] "Complete sequencing and characterization of 21,243 full-length human cDNAs." Ota T., Suzuki Y., Nishikawa T., Otsuki T., Sugiyama T., Irie R., Wakamatsu A., Hayashi K., Sato H., Nagai K., Kimura K., Makita H., Sekine M., Obayashi M., Nishi T., Shibaha

---

Note: This product is for in vitro research use only