

## 40S ribosomal protein S13 Polyclonal Antibody

Catalog No: #42356

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	40S ribosomal protein S13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 40S ribosomal protein S13 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 40S ribosomal protein S13 protein
Target Name	40S ribosomal protein S13
Other Names	40S ribosomal protein S13; ribosomal protein S13; RPS13; RS13
Accession No.	Swiss-Prot#: P62277
Uniprot	P62277
GeneID	6207;
Calculated MW	17kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting: □1:500 - 1:1000

Immunohistochemistry: 1:20 - 1:200

## Images



All lanes : 40S ribosomal protein S13 antibody at 2ug/ml

Lane 1 : EC109 whole cell lysate


Lane 2 : 293T whole cell lysate

Secondary

Goat polyclonal to Rabbit IgG at 1/15000 dilution

Predicted band size : 17kDa

Observed band size: 17kDa



Immunohistochemical analysis of paraffin-embedded human thyroid using #42356 at dilution of 1:100.

## Background

Belongs to the ribosomal protein S15P family.

## References

[1]"Cloning and analysis of the human S13 ribosomal protein cDNA." Chadeneau C., Lemoullac B., Denis M.G. Nucleic Acids Res. 21:2945-2945(1993) [2]"U14 snoRNAs are encoded in introns of human ribosomal protein S13 gene." Kenmochi N., Higa S., Yoshiham

Note: This product is for in vitro research use only