

General Concanavalin A (ConA) ELISA Kit

Catalog No: #EK13375



Package Size: #EK13375-1 48T #EK13375-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	General Concanavalin A (ConA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	General
Other Names	Con-A
Storage	12 months

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.065ng/mL

Sample Type:biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only