

StAR-related lipid transfer protein 7, mitochondrial Polyclonal Antibody

Catalog No: #42392

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	StAR-related lipid transfer protein 7, mitochondrial Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total StAR-related lipid transfer protein 7, mitochondrial polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human StAR-related lipid transfer protein 7, mitochondrial protein
Target Name	StAR-related lipid transfer protein 7, mitochondrial
Other Names	Gestational trophoblastic tumor protein 1 START domain-containing protein 7 STARD7 GTT1
Accession No.	Swiss-Prot#: Q9NQZ5
Uniprot	Q9NQZ5
GeneID	56910;
Calculated MW	43kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes: StAR-related lipid transfer protein 7, mitochondrial antibody at 2ug/ml + PC-3 whole cell lysate
 Secondary
 Goat polyclonal to Rabbit IgG at 1/10000 dilution
 Predicted band size: 43kDa
 Observed band size: 43kDa

References

[1] GTT1/StarD7, a novel phosphatidylcholine transfer protein-like highly expressed in gestational trophoblastic tumour: cloning and

Note: This product is for in vitro research use only