

# Protein phosphatase 1 regulatory subunit 11 Polyclonal Antibody

Catalog No: #42403

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Protein phosphatase 1 regulatory subunit 11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Protein phosphatase 1 regulatory subunit 11 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Protein phosphatase 1 regulatory subunit 11 protein
Target Name	Protein phosphatase 1 regulatory subunit 11
Other Names	PPP1R11
Accession No.	Swiss-Prot#: O60927
Uniprot	O60927
GeneID	6992;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human placenta using #42403 at dilution of 1:100.

## Background

Inhibitor of protein phosphatase 1.

## References

---

[1]"Cloning of a human homologue of the mouse Tctex-5 gene within the MHC class I region." Giffon T., Lepourcelet M., Pichon L., Jezequel P., Bouric P., Cam G., Pontarotti P., Gall J.-Y., David V. Immunogenetics 44:331-339(1996) [2]"Systematic sequenc

---

Note: This product is for in vitro research use only