

# Mouse D-Dimer (D2D) ELISA Kit

Catalog No: #EK13586



Package Size: #EK13586-1 48T #EK13586-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse D-Dimer (D2D) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mus musculus mouse
Other Names	D 2 Dimer
Storage	12 months

## Application Details

Detect Range:78.1-5000ng/mL

Sensitivity:25.4ng/mL

Sample Type:plasma.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only