

Rabbit D-Dimer (D2D) ELISA Kit

Catalog No: #EK13588



Package Size: #EK13588-1 48T #EK13588-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit D-Dimer (D2D) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Oryctolagus cuniculus Rabbit
Other Names	D 2 Dimer
Storage	12 months

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.27ng/mL

Sample Type:plasma.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only