

Transaldolase Polyclonal Antibody

Catalog No: #42412

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transaldolase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Transaldolase polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Transaldolase protein
Target Name	Transaldolase
Other Names	TALDO1
Accession No.	Swiss-Prot#: P37837
Uniprot	P37837
GeneID	6888;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human placenta using #42412 at dilution of 1:100.



Immunohistochemical analysis of paraffin-embedded human lung tissue using #42412 at dilution of 1:100.

Background

Transaldolase is important for the balance of metabolites in the pentose-phosphate pathway.

Note: This product is for in vitro research use only