

COMM domain-containing protein 4 Polyclonal Antibody

Catalog No: #42419

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COMM domain-containing protein 4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total COMM domain-containing protein 4 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human COMM domain-containing protein 4 protein
Target Name	COMM domain-containing protein 4
Other Names	COMMD4
Accession No.	Swiss-Prot#: Q9H0A8
Uniprot	Q9H0A8
GenID	54939;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human placenta using #42419 at dilution of 1:100.

Background

Down-regulates activation of NF-kappa-B. Ref.1

References

[1] "COMMD proteins, a novel family of structural and functional homologs of MURR1."Burstein E., Hoberg J.E., Wilkinson A.S., Rumble J.M., Csomos R.A., Komarck C.M., Maine G.N., Wilkinson J.C., Mayo M.W., Duckett C.S.J. Biol. Chem. 280:22222-22232(2005) [

Note: This product is for in vitro research use only